



## LDA AVATAR/FLOWSHEET

How to document line insertion, assessment and how to complete or discontinue an existing line


- To document line insertion go to the LDA Avatar on the top left of the flowsheets tab (1)



- (2) Click the +Add button and add your line, or select from (1)Mediport or (2) PIV



- Note: You can also add drains, wounds, etc...
- To document a implantable port and document it as a power port see below(3):



Enter a comment under port type by clicking the paper icon next to the port type (3)

Port Type

Single lumen port Double lumen port

4

Once you document a comment it will indicate there is text there by adding lines to the paper. (4)

**Implantable Port Right Chest**

Properties Orientation: Right Show all properties

Line Necessity

5

Once you have the established line you can click the “show all properties” hyperlink and will indicate this line has a comment associated with it. Hover to discover the comment “power port”

- Documenting assessment of a line (PIV or CVAD) under IV assessment (3), this includes unsuccessful insertion

Oncology Assessment IV Assessment Vital Signs Patient-Generated Sub... Patient Health Questi...

Accordion Expanded View All

1m 5m 10m 15m 30m 1h 2h 4h 8h 24h Interval Start: 0700 Reset Now

Appointment from 8/30/2023 i... 8/30/2023 0900

Search (Alt+Comma)

**Dominant Hand**

Which is your dominant hand? 6

**Peripheral IV 6/7/2023 Left;Proximal;Upper;Anterior Arm**

IV Properties Placed: 5/9/2023 0900 Show all

Site Assessment

Dressing Type

Line Status

Cap Change

Tubing Change

Dressing Status

Dressing Intervention

\$ Peripheral IV Charge

**Unsuccessful Insertion**

Name of Attempting Clinician

Number of Attempts

Line Type


Attempted Insertion Site (Vein)

Attempted Location

Additional Information

To document removal of a line **right click** the line header (6) Then select Remove

- One documentation is complete (removal was successful, date removed) it will show as below

|   |   |
|---|---|
| <b>[REMOVED] Peripheral IV 05/09/23 Left;Proximal;Upper;Anterior Arm</b> <span style="background-color: red; color: white; border-radius: 50%; padding: 2px 5px; font-weight: bold;">7</span> |   |
| IV Properties   |  <a href="#">Show all properties</a> |
| Site Assessment   |   |

- If you leave it as REMOVED it will remain on the flow sheet. To remove it from the flow sheet you must COMPLETE the IV. You COMPLETE the IV by **right clicking** the header a second time and selecting COMPLETE (7)