


# JOB AID – POCT DEVICES – Patient Identifier

## EPIC workflow change in patient identification from MRN to CSN for some POC devices

### CSN for the Patient Identifier:

Scan the square **2D barcode** found on the patient wristband with these devices:

- Nova StatStrip:            Glucose            All locations
- TEG 6s:                    TEG                OR locations only
- GEM 5000:                Blood gas            OR locations or emergent patients only
- CoaguChek                PT/INR            Outpatient locations only – AMB patient label
- DCA Vantage              HbA1c             Outpatient locations only – AMB patient label
- Afinion                    HbA1c             Outpatient locations only - AMB patient label
- Clinitek                    UA                 Outpatient locations only – ABM patient label



### MRN for the Patient Identifier:

Scan the **linear MRN barcode** found on the patient wristband or patient label with these devices:

- Istat:                      Creatinine
- Nova StatSensor:        Creatinine
- Signature Elite:         ACT

