


JOB AID – POCT DEVICES – Patient Identifier

EPIC workflow change in patient identification from MRN to CSN for some POC devices

CSN for the Patient Identifier:

Scan the square **2D barcode** found on the patient wristband with these devices:

- Nova StatStrip: Glucose All locations
- TEG 6s: TEG OR locations only
- GEM 5000: Blood gas OR locations or emergent patients only
- CoaguChek PT/INR Outpatient locations only – AMB patient label
- DCA Vantage HbA1c Outpatient locations only – AMB patient label
- Afinion HbA1c Outpatient locations only - AMB patient label
- Clinitek UA Outpatient locations only – ABM patient label



MRN for the Patient Identifier:

Scan the **linear MRN barcode** found on the patient wristband or patient label with these devices:

- Istat: Creatinine
- Nova StatSensor: Creatinine
- Signature Elite: ACT

